



Health Care Call System

Geriatric Callcord

Cat. 6879

DESCRIPTION

Jeron Geriatric Callcords are used in Jeron Nurse Call Systems to provide a means of signaling for assistance with only slight hand or body pressure. Geriatric Callcords are used in hospitals and nursing homes for geriatric or pediatric patients, for feeble or disabled patients, for complete safety in an oxygen-enriched environment, or for use in showers.

OPERATION

Slight hand or body pressure on the bulb places a call on the system as described on the system and patient station data sheets.

Removal of the callcord from the patient station produces a cordout call which can only be canceled by replacing the cord into the jack. Cat. 6836 Dummy Plug may be inserted into the jack instead of the callcord to preclude a cordout call.



EQUIPMENT

Geriatric Callcords have a gray, non-toxic tubing with matching molded bulb for use in oxygen enriched environments. The cord terminates in a molded plastic plug with an 8-pin connector. A plated spring-type metal bed clip attaches the callcord to the bed. Bulb, cord, and plug withstand common hospital cleaning agents.

TECHNICAL DATA

Bulb and Cord:

- 6 ft. (1.8 m) gray non-toxic tubing with matching molded air bulb

Controls:

- Air Bulb; pressure activated

Connector:

- 9-pin plug; molded plastic shell

Certification:

- UL 1069 Listed
- C-UL — Canadian Standard C22.2 NO.125

PRODUCT IS SUBJECT TO CHANGE WITHOUT NOTIFICATION.