



Health Care Call System

Terminal Block

Cat. 6171

DESCRIPTION

Jeron's Cat. 6171 Terminal Block is used in Pro-Alert™ 610 Tone-Visual Nurse Call Systems to interconnect the Cat. 6150/6153 Main Control with room station field wiring.

EQUIPMENT

The Cat. 6171 Terminal Block contains two rows of two-clip terminals which are pre-wired to two 50 pin (25-pair) female connectors, side mounted for top cable entry. Each Cat. 6171 block has the capacity for 96 room stations. The Terminal Blocks are furnished with two, ten-foot, 25-pair plug-ended cables for connection to the Cat. 6150/6153 Main Control.

For 48 Room Station capacity use the Cat. 6671 Terminal Block



TECHNICAL DATA

Mounting:

- Built-in mounting brackets for surface mounting

Dimensions (H x W x D):

- 10 x 3 x 2.75 in.
(24.4 x 7.6 x 7 cm)

Capacity:

- 96 Room Stations

Field Connection:

- Requires Cat. 8880 Punch Block Tool for wiring terminal clips

Backbox (flush mount only):

- Cat. 8888, 15 x 12 x 4 in. screw cover

Certification:

- UL 1069 Listed
- C-UL—Canadian Standard C22.2 NO.125

PRODUCT SUBJECT TO CHANGE WITHOUT NOTIFICATION.